



## **NEW** VS series steam sterilizers

A full range of compact, efficient, configurable autoclaves focused on process safety, with the widest choice of loading options!

Featuring the

## **4D<sup>IR</sup> SENSOR**

For state of the art direct monitoring of steam **saturation** & **penetration** of every sterilization load!

Find more details:

### **Stand no. 5**

ground floor

### **Satellite Symposium**

Yangtze 2 room

Friday, November 1st, 12:30 AM - 2:00 PM

[www.steelcogroup.com](http://www.steelcogroup.com)

**Miele**

Group  
Member

**NEW**

# VS series - futureproof sterilizers

## Advanced Technology

Efficiency  
& Sustainability



**ECO OPTIONS**

minimized  
water & energy consumptions

Process Safety

**4D**

**IR SENSOR**

direct monitoring  
of steam saturation & penetration

**NCG SYSTEM**

measurement  
of Non Condensable Gasses

Performance  
& Reliability



**DESCALING**

automatic descaling  
of the vacuum pump

## Flexible Configuration



from 4 to 12 STUs  
**compact footprint**  
950 mm width, full frontal access



from 4 to 12 STUs  
**right or left side**  
**HMI and technical area**



from 6 to 18 STUs  
suitable **for logistic systems**

## Loading Automation



Automatic loading  
transfer trolleys



ATS - fully automated  
transfer loading system



Easy Load system  
for stackable loads